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**Definition**

Onlays/Occlusions are polymerized, laser-treated, prefabricated, radio-paque, highly filled nano-hybrid composite enamel occlusion shells. The occlusal and anatomical profile templates in form of a thin composite shell facilitate the direct veneering of individual or more posterior teeth. Onlays/Occlusions are quick and easy and provides an esthetic, functional and high quality outcome.

Onlays/Occlusions are cemented with the same material they are made of, which leads to a stable unit of the entire restoration (between dentine/enamel + composite + onlays/occlusions).

Onlays/Occlusions are made of a single and complex posterior tooth restorations with a biomechanical result. The pre-fabricated Onlays/Occlusions are the anatomical basis for single or complete reconstructions as well as for lifting the vertical dimension in the posterior region. This also makes it possible to attain functioning guidance of anterior teeth by using edelweiss DIRECT VENEARS.

**Composition:**

Barium dental glass
Bis-GMA-based matrix made of dental resins
Pigments additives catalyst

**Technical data:**

Filler content: 82 wt %
Variation range of inorganic filler particles lies between 0.02-3 µm.

**Features of edelweiss COMPOSITE**

(which is recommended for mounting and direct filling therapy)

- Excellent handling
- Excellent processing properties
- Easy and fast polishing
- Low shrinkage due to nano-technology and high degree of filling (82%)
- Excellent abrasion behavior
- After fixation: antibacterial properties
- Zinc and fluorine component provides antibacterial surface in filler
- Natural fluorescence

**Indications:**

- Caries
- Old restorations
- All filling classes 1-4
- Misaligned teeth
- Tooth fractures
- Tooth discoloration
- Discoloration
- Veneering of old crowns
- Anatomical malformation
- Abrasion
- Erosion
- Cosmetic corrections
- Lifting of vertical dimension in the posterior region

**Recommended application**

- Clean the tooth to be restored and the adjacent mesial and distal teeth with fluoride-free cleaning paste. Interdental roughing/cleaning with abrasive strips.
  - Selection of tooth shade
    - The tooth shade is selected using the supplied sizing guide (edelweiss OCLUSION ↗ Sizing Guide). The sizing guide is placed over the teeth to be restored, whereby the contours visible on the transparent film are used to select the best matching edelweiss OCLUSION ↗.
- Selecting the color shade
  - Determine the dentine shade in the buccal & occlusal region with the dentine shade tip as the shade is most pronounced here. Place dentine tip in the Enamel VENEER tray. (first add glycerine to the tray for light refraction.)
- Selection of shade brightness
  - Select shade of composite, in case of severely damaged teeth or discolored base, a dentine material is recommended to achieve homogeneous coloring. Enamel composite can be used to correct the shade Enamel.
- General pretreatment
  - Drying (coffer dam)
  - Working under dry conditions is the basis for optimal results.
- As minimally invasive as possible, depending on indication.

- Matrices/interdental wedges
  - Apply a thin, transparent matrix band to protect the adjacent teeth and for separation purposes. If required, fixate the matrix approxi-mally with interdental wedges.
- Pretreatment edelweiss OCLUSION ↗
  - edelweiss OCLUSION ↗ are fitted, prepared accordingly, and cleaned with phosphoric acid etching gel depending on the indication.
- 2) Improve bonding, the back of the edelweiss OCLUSION ↗ is sand-blasted with 25µm or 50µm aluminum oxide or roughened with a diamond composite instrument.
- 5.1 edelweiss OCLUSION ↗ are fitted, prepared accordingly, and cleaned with phosphoric acid etch gel and allow to work for at least 30s to 60s on enamel, resp. 15s on dentine.
- 3) Rinse thoroughly for approx. 20s and air blow gently with oil and water-free compressed air.
7. Mounting of edelweiss OCLUSION ↗ (direct application)
  - The use of edelweiss nano-hybrid composite is recommended to ensure esthetic and functional success.
- Dispense Dentin/Enamel Bond (Peak univ., Optibond, etc.) directly on the Microbrush and apply it to the occlusal region free-hand.
- 7.2 Air blow briefly, Preharden with light output of less than 1000mW/cm² for 30s.
- 7.3 Place matrix bands, (i.e. Triodent, Hawe), if required, build up the cavity resp. lingually and buccally with several layers of edelweiss Composite and instead of modeling the occlusal region free-hand, composite is applied to the edelweiss OCLUSION ↗ (s bonded) and positioned on the teeth to be treated, excess is removed and fixated in final position with a UV lamp. Remove matrix bands. Remove excess with rotating instrument. Avoid adhesive puddles as these may impair the accuracy of fit of the edelweiss OCLUSION ↗.
- 5.7 Then light-cure the edelweiss VENEER Bond for at least 20s.
6. Pretreatment of teeth
  - The use of a coffer dam is recommended.
- 6.1 Clean tooth preparation
- 6.2 Etch with 37% phosphoric acid etch gel and allow to work for at least 30s to 60s on enamel, resp. 15s on dentine.
- 6.3 Rinse thoroughly for approx. 20s and air blow gently with oil and water-free compressed air.
7. Mounting of edelweiss OCLUSION ↗ (direct application)
  - The use of edelweiss nano-hybrid composite is recommended to ensure esthetic and functional success.
- Dispense Dentin/Enamel Bond (Peak univ., Optibond, etc.) directly on the Microbrush and apply it to the occlusal region free-hand.
- 7.2 Air blow briefly, Preharden with light output of less than 1000mW/cm² for 30s.

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**Indirect application:**

By lining the pre-fabricated occlusal surfaces on existing model occlusal with the nano-hybrid composite, one can simply create a new functioning biomechanical, balanced and smoothly working mastilcular apparatus without irritations. As a consequence, the occlusions of the model are transferred to the patient’s occlusion, fixed with bonding and cemented with composite.

**Direct procedure, but for advanced professionals for one-sided (mandibular) lifting of the vertical dimension:**

- The maxilla and the mandible models are articulated facially-adequate in habitual teething. Then protrusion is performed as required and the supporting pin of the articulator raised. Partial bite splints with “stops” are prepared for the mandibular quadrants. These are fixed temporarily in the mouth, in other words, provisionally. Already after a few days this will represent the ideal position of the mandible, possibly requiring reocclusion. On the basis of these bite splints the mandible is articulated anew versus the maxilla.

Instead of lining the occlusions on the model in the articulator, a bite registry material (silicone) is used to make an occlusal impression which accurately portrays the clearance between maxilla and mandible in the articulator. Using the fixed silicone bite register, direct application of the pre-fabricated occlusion (in the patient’s mouth) can be performed. The silicone bite register is in each case placed unilaterally (left or right) on the patient’s occlusion. In other words, when the left quadrant is restored with occlusions, the right quadrant fixates the height and position of bite elevation and thus the final position of the composite-lined and bonded occlusions on the left when biting down using the silicone bite registry (replaces the supporting pin when closing the articulator).

Once the left quadrant has been restored, the right quadrant is treated. After fixing the edelweiss VENEER Bond in the parte posterior/parte interior basal processeda of edelweiss OCLUSION ↗. Aplicar una gota de edelweiss VENEER Bond en la edelweiss OCLUSION ↗ con un pincel y extender durante al menos 30s.

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- None known if used as intended.

**Side effects**

Systemic side effects are not known. Individual cases of contact allergies have been reported for products with a similar composition in cured condition.

**Interactions**

Discoloration may occur during contact with cationic mouth washes as well as plaque releators or chlorhexidine.

**Shelf-life and labeling**

The expiry date and the batch description are given on the container. Do not use after expiry date.

**Note:**

May only be supplied to dentists and dental laboratories or on their behalf. Keep away from children! Only for use by dentists!

**Important**

To ensure an optimal result, it is recommended to use in combination with all edelweiss components.

Date of issue of this package information leaflet 2018-03

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