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Instruction for use English

Definition

edelweiss Nano-Hybrid Composite is a light curing, radiopaque, highly filled nano-hybrid composite filling material for restorations that meet the highest standards.

The system consists of 2 components:

Dentin in 5 shades vita shade (A0 to A3.5)
Enamel

edelweiss Nano-Hybrid Composite is light cured directly in the tooth cavity or on a model using any standard curing light (e.g. halogen, LED).

Composition

Barium dental glass
Bis-GMA based resin matrix
Pigments
Additives
Catalyst

Technical data

Filler ratio: 82 % by weight = 65 % by volume. The variation of inorganic filler particle is between 0.02-3 µm.

Indications

Direct applications:
– Cementation of edelweiss VENEER
– Filling of cavity classes I, II, III, IV, V

Indirect applications:

- Inlays
- Onlays
- Veneers

Dentin is specially indicated for:
– Reconstruction of natural dentin

Enamel is specially indicated for:
– Reconstruction of natural enamel

Recommended application

- Clean the tooth surface prior to preparation. Remove any cleaning residue using water spray. Select the shade using the edelweiss Dentin & Enamel Guide. The use of a dental rubber dam is recommended.
- Prepare the cavity according to the general requirements for adhesive dentistry. Bevel the enamel margins. Protect any areas near the pulp.
- If necessary, etch according to your preferred etching technique, e.g. Dentin/Enamel Bond (Peak Univ., Optibond, etc.). Follow the instructions of the material used.
- Apply adhesive according to its instructions for use.
- Apply a liner using a flowable composite or compomer, as necessary. The use of a matrix is recommended.
- Apply the desired shade of the edelweiss Nano-Hybrid Composite into the cavity using a maximum thickness layer of 2 mm. In all cases, the recommended curing time for each layer is 20 s.
- Remove excess, contour and finish, using e.g. diamond finishing burs and flexible discs in decreasing grit sizes. A high gloss finish for the edelweiss Nano-Hybrid Composite can be achieved using polishing paste together with polishing cups/discs.

Essential notes

- Light curing units should have an output of 450 nm and should be checked regularly. The curing light output intensity should be at least 400 mW/cm². Place the curing light as close as possible to the restorative material.
- Do not remove the oxygen-inhibition layer, since this is needed to create a sufficient bond to the next layer.
- Follow the instructions for use for all materials used together with edelweiss Nano-Hybrid Composite.

Contraindications

Application of edelweiss Nano-Hybrid Composite is contraindicated if a dry working area or the recommended application technique is not possible, if the pulp is exposed, or if there are any known allergies to any components in the edelweiss Nano-Hybrid Composite.

Side effects

Side effects are not known to date. In individual cases, hypersensitivity to components of the material cannot be excluded. Avoid contact of the paste with skin. If accidental contact occurs, wash immediately using soap and water. In case of contact with eyes, rinse immediately using plenty of water, and consult a physician.

Interactions

Avoid materials containing eugenol, moisture or oily air, since they can inhibit polymerization at the contact area. The use of cationic mouth-washes, plaque indicators or chlorohexidine may lead to discolorations.

Classification

edelweiss Nano-Hybrid Composite fulfills the requirements of ISO 4049:2000. The specified curing depth of 2 mm (20 s) for edelweiss Nano-Hybrid Composite is prepared according to ISO 4049:2000. The specified thickness layer of 2 mm and light curing time of 20 s for each layer is recommended for all cases.

Storage

Store at normal room temperature (15-25 °C / 59-77 °F). Please close syringe / composite tip immediately after each use to prevent premature setting from ambient light.

Shelf life and marking

Expiration date and the lot number are marked on the immediate containers and packaging. Do not use after the date of expiration.

Note

Only to be supplied to dentists and dental laboratories or upon their request. Keep out of the reach of children! For dental use only!

Important

To ensure for optimal clinical results, we recommend that all edelweiss components in this system are used in conjunction with each other.

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