


INDICATION GUIDE

DISCOVER YOUR SMILE

edelweiss
DENTISTRY

shaping the future of dentistry



BRINGING HARMONY AND BALANCE
TO THE FACE NATURALLY

BIO-ESTHETICS AND FUNCTION IN ONE APPOINTMENT

edelweiss CAD/CAM BLOCK T-BLOCK/C-BLOCK/i-BLOCK

The new edelweiss CAD/CAM BLOCKs have a homogenous glass-ceramic phase manufactured through a process of controlled laser sintering technology fusing the individual glass crystals. For the first time, the properties of both particle reinforced and predominately glass-ceramics are featured in a single hybrid glass block, having the strength of reinforced ceramics, the optical properties of predominantly glass ceramics and the flexibility of dentin. Its versatile application together with its time and cost saving procedure makes the edelweiss CAD/CAM BLOCKs a sound investment in the future with the patient's best interest in mind.

Convince yourself

THE EDELWEISS DIGITAL WORKFLOW

T-BLOCK

Translucent Enamel



i-BLOCK Implant

Translucent Enamel & Chroma



C-BLOCK

Chroma



A0



A1



A2



A3

The new era in CAD/CAM design by edelweiss dentistry, combines modern laser sintering to produce a solid hybrid glass phase to achieve optimum esthetics and high strength.

Three technologies in a single edelweiss CAD/CAM BLOCK:

- glass-ceramics for esthetics
- particle filled ceramics for strength
- resin technology for resilience



The difference lies in it's manufacturing process.

ATTRITION / ABRASION

edelweiss CAD/CAM BLOCK





MISALIGNED / DIASTEMA

edelweiss CAD/CAM BLOCK






DISCOLOURED / CHIPPED TEETH

edelweiss CAD/CAM BLOCK





DISCOLOURED / SHORTENED TEETH

edelweiss CAD/CAM BLOCK





DISCOLOURED / CHIPPED TEETH

edelweiss CAD/CAM BLOCK



CROWN / VENEERS

edelweiss CAD/CAM BLOCK



DISCOLOURED TEETH / DIASTEMA

edelweiss CAD/CAM BLOCK



TOOTH WEAR / DISCOLOURED TEETH

edelweiss CAD/CAM BLOCK



ATTRITION / ABRASION

edelweiss CAD/CAM BLOCK





BRUXISM / TOOTH WEAR

edelweiss CAD/CAM BLOCK





RETROCLINED TEETH

edelweiss CAD/CAM BLOCK





RETROCLINED TEETH

edelweiss CAD/CAM BLOCK





RETROCLINED TEETH

edelweiss CAD/CAM BLOCK



BRUXISM / TOOTH WEAR

edelweiss CAD/CAM BLOCK



DISCOLOURED TEETH

edelweiss CAD/CAM BLOCK





i-BLOCK: Customised screw-retained implant supported crowns



edelweiss CAD/CAM BLOCK



*Single unit as
abutment crown*




*Mesostructure
and crown*



i-BLOCK



edelweiss CAD/CAM BLOCK





i-BLOCK – the procedure: Model based digital workflow


edelweiss CAD/CAM BLOCK



Surgical placement of implant and intra oral scan of the initial situation.

The base edelweiss iBLOCK + crown (edelweiss C-BLOCK) was designed and milled.
The TiBase + edelweiss Base + edelweiss Crown was checked for fit.
Final placement and checked in occlusion.

Immediate loading digital workflow



Initial situation scanned and virtual extraction done on patient.
Implant and crown planned virtually.
On the patient, tooth was extracted and implant placed.

Custom base was milled with edelweiss iBLOCK for optimum biological seal and emergence profile.
The crown was milled with edelweiss C-BLOCK.
The TiBase + edelweiss Base + edelweiss Crown was fit together and placed onto the implant.

Step by Step

edelweiss CAD/CAM BLOCK

TOOTH PREPARATION



1

Pre-op, non-prep.



2

Try in of veneer.



3

Etching with 37 % phosphoric acid.




4

Rinse thoroughly with water to remove etchant.



5

Application of adhesive bonding agent to tooth surface.



6

Light-cure for 20 seconds.



7

Placing a thin layer of edelweiss Bond on the veneer and light-cure for 10 seconds.



8

Evenly adapt the Nano-Hybrid Composite within the veneer using an instrument (spatula).




9

Placing veneer.



10

Removal of excess material at the marginal area.



11

Light-cure for 20 seconds.



12

Finishing



13

Retraction cord removal.



14

Finishing strip



15

Post-op



shaping the future of dentistry

edelveiss dentistry ® and the edelveiss Logo ® are registered trademarks of edelveiss dentistry products gmbh • Austria

edelveiss dentistry products gmbh • Austria
office@edelveissdentistry.com
www.edelveissdentistry.com



822300 / 2026-03



© 2026 edelveiss dentistry products gmbh • Austria. All material contained in this brochure are the property of edelveiss dentistry products gmbh • Austria. No part of this publication may be reproduced, distributed, or transmitted in any form or by any means, without the prior written permission of the edelveiss dentistry products gmbh • Austria. All Rights Reserved.